

AUTOMATIC DENSITY METER CODE MDY-MD30

DATA
OUTPUT

- According to GB/T 29617, GB/T 2013, ISO 12185, ASTM D4052, ASTM D5002, ASTM D3505
- Rapid measurement of liquids such as sugar solutions, fruit juices, alcohols, oils, food additives, metalworking cleaners, etc.
- 2~3 mL of sample is needed each time
- Built-in Peltier temperature control device for more accurate measurements
- Dry air and distilled water can be used for daily calibration
- Compatible with automatic sampling and manual sampling
- Rapid pre-test mode and precision mode
- Data can be saved and reports can be exported



manual sampling



automatic sampling



drying barrel (included)

SPECIFICATION

| | |
|-------------------------------|---|
| Density range | 0.0000g/cm ³ ~3.0000g/cm ³ |
| Measurement method | density g/cm ³ , density kg/m ³ , specific gravity of water |
| Density resolution | 0.0001g/cm ³ |
| Density accuracy | ±0.0005g/cm ³ |
| Density repeatability | ≤0.0002g/cm ³ |
| Temperature resolution | 0.01°C |
| Temperature control range | 5°C~50°C |
| Temperature control precision | ±0.05°C |
| Temperature control accuracy | ±0.1°C |
| Observation | high definition camera/video observation |
| Data storage | 128G |
| Display | 10.1 inch color LCD touch screen |
| Power supply | 100V~240V, 50Hz~60Hz |
| Dimension | 439×350×217mm |
| Weight | 9kg |

STANDARD DELIVERY

| | |
|----------------------------|------|
| Main unit | 1pc |
| Power adapter | 1pc |
| PP syringe (MDY-MD30-ZSQ1) | 3pcs |
| Drying barrel | 1pc |
| Waste liquid bottle | 1pc |
| Silicone tube connector | 1set |
| Silicone tube | 1pc |
| Dust cover | 1pc |

OPTIONAL ACCESSORY

| | |
|-----------------|--------------------|
| Stylus printer | MDY-MD30-PRINTER-A |
| Thermal printer | MDY-MD30-PRINTER-B |
| PTFE syringe | MDY-MD30-ZSQ2 |
| Glass syringe | MDY-MD30-ZSQ3 |